


FEHLING...
... THE DIFFERENCE

8

ASPIRATORS
SAUGER



FEHLING
SURGICAL INSTRUMENTS, INC.

509 Broadstone Lane · Acworth, GA 30101


Toll Free: 1-800-FEHLING (800-3345464), Fax: 770 794-0122, E-Mail: info@fehlingsurgical.com



MALLEABLE / BIEGBAR

Ø TUBE / ROHR	LENGTH / LÄNGE A						
	80 mm	100 mm	125 mm	145 mm	165 mm	180 mm	220 mm
0,6 mm	HAU-0	-	-	-	-	-	-
0,8 mm	HAU-4	-	-	HAI-0	-	-	-
1,0 mm	HAU-5	HAU-1	HAU-6	HAI-7	HAX-4	-	-
1,2 mm	-	HAU-2	HAU-7	HAI-8	-	-	-
1,6 mm	-	HAU-3	HAU-8	HAI-9	-	-	-
1,8 mm	-	-	-	-	HAI-1	HAI-4	HAX-1
2,0 mm	-	-	-	HAU-9	HAI-2	HAI-5	HAX-2
2,3 mm	-	-	-	-	HAI-3	HAI-6	HAX-3
3,0 mm	-	-	-	-	HAX-5	-	-

MICRO ASPIRATORS
MIKROSAUGER



MALLEABLE / BIEGBAR

Ø TUBE / ROHR	LENGTH / LÄNGE A	
	120 mm	160 mm
1,0 mm	VOA-5	-
2,0 mm	VOA-6	VOA-1
3,0 mm	VOA-7	VOA-2
4,0 mm	VOA-8	VOA-3
5,0 mm	VOA-9	VOA-4
6,0 mm	-	VOA-0

CERAMO® SUCTION CANULAS
CERAMO® SAUGKANÜLEN

FRAZIER / FERGUSSON



CERAMO®

Ø TUBE / ROHR	TOTAL LENGTH / GESAMTLÄNGE	
	180 mm	200 mm
1,0 mm	–	VOC-1
2,0 mm	VOB-2	VOC-2
2,5 mm	VOB-3	VOC-3
3,0 mm	VOB-4	VOC-4
4,0 mm	VOB-5	VOC-5
5,0 mm	VOB-6	VOC-6



UNCOATED
UNBESCHICHTET

Ø TUBE / ROHR	TOTAL LENGTH / GESAMTLÄNGE		
	180 mm	200 mm	250 mm
1,0 mm	–	HAR-1	–
2,0 mm	NAT-7	HAR-2	NHR-8
2,5 mm	NAT-8	HAR-0	NHS-0
3,0 mm	NAT-9	HAR-3	NHR-9
4,0 mm	NAU-1	HAR-4	NHS-1
5,0 mm	NAU-2	–	NHS-2

Ø 3,5 / 1,5 mm



JAMIESON
DISSECTING ASPIRATOR
DISSEKTIONSSAUGER
270 mm
MMA-7

STRAINER ONLY
SAUGKORB ALLEIN
MHE-7



HLM
ASPIRATOR
SAUGER
280 mm
MEL-7

DECALCIFICATION APSIRATOR
ENTKALKUNGSSAUGER



Ø 6,0 mm

$\frac{2}{3}$

240 mm

MNF-5

ASPIRATORS FOR PEDIATRIC SURGERY
SAUGER FÜR KINDERCHIRURGIE

Ø 3,3 mm



210 mm
MNE-5

Ø 4 mm



225 mm
MCE-9

Ø 4 mm



245 mm
MCE-8

GRECO

ASPIRATORS FOR PEDIATRIC SURGERY
SAUGER FÜR KINDERCHIRURGIE



Ø TUBE / ROHR	210 mm
3,3 mm	MNL-2
5,0 mm	MNL-5

Ø TUBE / ROHR	180 mm
3,3 mm	MNL-1
5,0 mm	MNL-4

Ø TUBE / ROHR	150 mm
3,3 mm	MNL-3
5,0 mm	MNL-6

GRECO

ASPIRATORS FOR PEDIATRIC SURGERY
SAUGER FÜR KINDERCHIRURGIE



DETACHABLE BASKET / KORB ABNEHMBAR

Ø TUBE / ROHR	180 mm
3,3 mm	MNM-1
5,0 mm	MNM-3



DETACHABLE BASKET / KORB ABNEHMBAR

Ø TUBE / ROHR	160 mm
3,3 mm	MNM-2
5,0 mm	MNM-4

GRECO

ASPIRATORS FOR PEDIATRIC SURGERY
SAUGER FÜR KINDERCHIRURGIE



DETACHABLE BASKET / KORB ABNEHMBAR


Ø TUBE / ROHR	230 mm
6,0 mm	MNL-7

DETACHABLE BASKET / KORB ABNEHMBAR


Ø TUBE / ROHR	205 mm
6,0 mm	MNL-8

GRECO

HLM-CANNULAS FOR PEDIATRIC SURGERY
HLM-KANÜLEN FÜR KINDERCHIRURGIE



∅ TUBE / ROHR	95 mm
3,3 mm	MNM-5
5,0 mm	MNM-7



∅ TUBE / ROHR	80 mm
3,3 mm	MNM-6
5,0 mm	MNM-8



Ø 7 mm
L 25 mm
MKK-5

Ø 6 mm
L 23 mm
MKN-8

PEDIATRIC VENTRICLE ASPIRATOR
PÄDIATRISCHER VENTRIKELSAUGER



ICE WATER ASPIRATOR
EISWASSERSAUGER
MCC-1



Ø 2,5 mm
L 48 mm
HAH-9

Ø 3,7 mm
L 55 mm
HAJ-2

Ø 5,0 mm
L 160 mm
HAJ-3

HEPARIN CANNULA
HEPARINKANÜLE



AORTIC DEAIRING CANNULA
AORTEN-ENTLÜFTUNGSKANÜLE
MBR-2

BUTTONED CANNULAS
KNOPFKANÜLEN

MALLEABLE
BIEGBAR



Ø 1,2 mm
LENGTH / LÄNGE
80 mm
VAD-3



Ø 1,5 mm
LENGTH / LÄNGE
80 mm
VAD-5



Ø 1,2 mm
LENGTH / LÄNGE
80 mm
VAD-4



Ø 1,5 mm
LENGTH / LÄNGE
80 mm
VAD-6